# **Endless Drive Progressive**



# **Endless Drive Progressive**

### A premium quality progressive lens for day or night driving

A standard progressive lens may be uncomfortable to driving. Perhaps the visual field is insufficient, or peripheral vision is limited. At night, the lack of light complicates matters, increasing hazards and causing fatigue. In addition, one out of three wearer experience night myopia, considerably reducing their ability to focus on distant objects.

Endless Drive Progressive is a unique customized lens. Its power distribution is specifically adapted for driving. Its night vision zone compensates for the difference in refractive error that occurs between day and night. It provides greater visual acuity and reduces eyestrain.

By incorporating IOT Digital Ray-Path® 2, Technology the lens optimization process takes into account the wearer's accommodative ability to focus on different distances. The level of customization is therefore much more sophisticated and effective than previous methods.

97% of wearers who have used **Endless Drive Progressive** lenses have felt much more comfortable driving both during the day and at night.

Endless Drive Progressive lenses can be recommended for general use.

### Category

Use	Special
Product	Personalized
Frequency of use	Regular

### **Ideal** wearer

Drivers looking for the most innovative solutions.

Wearers looking for greater comfort, visual quality, and safety while driving.

Wearers who feel they see worse at night than during the day.

Wearers with all types of **prescription and** addition powers.




### Unique power distribution based on the visual needs of drivers

Endless Drive Progressive lenses have a unique power distribution based on the way we see while driving. They have been developed to maximize the near and mid-range vision zones.

Its unique power distribution is based on three key aspects:

#### **Endless Drive Progressive**

- Has a night vision area.
- · The distance area is larger.
- · The higher area is wider.
- · Lateral astigmatism is considerably reduced.



Progression profile: This lens incorporates extra negative power (up to 0.25 D) in the upper area of the lens to compensate for the effects of night myopia.

Unwanted astigmatism: Endless Drive Progressive lenses have a maximum astigmatism level 14.6% lower than that of general use premium progressive lenses. This improves vision in the peripheral areas of the lens and increases the visual range.

Fitting Cross: Endless Drive
Progressive lenses are 70% wider in
the distance and 45% in the upperintermediate when compared to
premium general-use lenses. Wearers
will enjoy expanded views of the

road, dashboard, and mirrors.



### Infinite design configuration options

Create a unique product suited to your market needs

**Endless Drive Progressive** lenses allow for countless configurations, providing the opportunity to differentiate for unique market needs.

No other lens design portfolio is as modular, flexible, and versatile as IOT's.

VALUES	POSSIBILITIES		
MFHS & INSETS			
Minimum fitting height (MFH)	18 mm		
Insets	Automatic or manual (from 0 to 4 mm – infinitesimal steps)		
POWER RANGES		Compatibili	$\mapsto$
Rx range	Extended to your blank limits		
Add powers	From 1.00 D to 4.50 D (infinitesimal steps)	Material & blank	Endless Drive Progressive lenses are compatible with any blank provider and lens index.
THICKNESS		provider	
Thinning prism	Standard or equalized	Castings	Endless Drive Progressive lenses are compatible with any coatings you run at your lab.
Thickness optimization	Decentration or lenticularization	Coatings	
LENS COMPENSATION			
Personalization parameters	Real or defaults	Machinery & LMS	Endless Drive Progressive lenses are compatible with almost any machinery supplier and LMS.
Power compensation	Optimal or customized		
Prism compensation	Enable or disable		Supplier and EMS.
REFERENCE POINTS			
Layout	Standard or on demand		
MANUFACTURING			
Crib	Rounder/elliptical/shaped		
Prism	Blocked, generated or mixed		

## **Endless Drive Progressive**

### **Technologies**

DRP 2

### Features —



#### Personalization

The back surface is Personalized according to the use parameters, creating a unique lens for each wearer. If no actual parameters are available, IOT Digital Ray-Path® 2 uses them by default.



#### Compensated power

Lens power differs from prescribed power. The design is calculated, point by point, to ensure wearers perceive the proper power when looking through their lenses at every distance and direction of gaze, including near distances for viewing electronic devices.



# Optimized in accommodative space

IOT Digital Ray-Path® 2 incorporates the intelligent use of the wearer's accommodation into the traditional calculations for reducing oblique aberrations, resulting in a superior personalized lens. Oblique aberrations are minimized more effectively than ever before.



### Consistency

The perceived power distribution remains stable, regardless of the prescription or base curve. This is especially beneficial for high prescriptions and large or wrapped frames.

